

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revK_14Aug2019
 Last update: 27/12/2019



Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 500 ng polyA-tailed RNA in 9 µl	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
<input type="checkbox"/> Direct RNA Sequencing Kit (SQK-RNA002)	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937)	<input type="checkbox"/> Microfuge
	<input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> SuperScript III Reverse Transcriptase (ThermoFisher Scientific, 18080044)	<input type="checkbox"/> Ice bucket with ice
	<input type="checkbox"/> 10 mM dNTP solution (e.g. NEB N0447)	<input type="checkbox"/> Timer
	<input type="checkbox"/> NEBNext Quick Ligation Module (E6056)	<input type="checkbox"/> Thermal cycler
	<input type="checkbox"/> Agencourt RNAClean XP beads	<input type="checkbox"/> Qubit fluorometer (or equivalent for QC check)
	<input type="checkbox"/> Qubit RNA HS Assay Kit (ThermoFisher Q32852)	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> Qubit dsDNA HS Assay Kit (ThermoFisher Q32851)	

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Library preparation</p> <p>Prepare the RNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 500 ng RNA into a 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Adjust the volume to 9 µl with Nuclease-free water <input type="checkbox"/> Mix thoroughly by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <p>In a 0.2 ml thin-walled PCR tube, mix the reagents in the following order:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 3.0 µl NEBNext Quick Ligation Reaction Buffer <input type="checkbox"/> 9.0 µl RNA <input type="checkbox"/> 0.5 µl RNA CS (RCS), 110 nM <input type="checkbox"/> 1.0 µl RT Adapter (RTA) <input type="checkbox"/> 1.5 µl T4 DNA Ligase <ul style="list-style-type: none"> <input type="checkbox"/> Mix by pipetting and spin down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. 	

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revK_14Aug2019
 Last update: 27/12/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Mix the following reagents together to make the reverse transcription master mix:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 9.0 µl Nuclease-free water <input type="checkbox"/> 2.0 µl 10 mM dNTPs <input type="checkbox"/> 8.0 µl 5x first-strand buffer <input type="checkbox"/> 4.0 µl 0.1 M DTT <p><input type="checkbox"/> Add the master mix to the 0.2 ml PCR tube containing the RT adapter-ligated RNA from the "RT Adapter ligation" step. Mix by pipetting.</p> <p><input type="checkbox"/> Add 2 µl of SuperScript III reverse transcriptase to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Place the tube in a thermal cycler and incubate at 50° C for 50 min, then 70° C for 10 min, and bring the sample to 4° C before proceeding to the next step.</p> <p><input type="checkbox"/> Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Resuspend the stock of Agencourt RNAClean XP beads by vortexing.</p> <p><input type="checkbox"/> Add 72 µl of resuspended RNAClean XP beads to the reverse transcription reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 200 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p>Keep the tube on magnet, and wash the beads with 150 µl of freshly prepared 70% ethanol without disturbing the pellet as described below.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Keeping the magnetic rack on the benchtop, rotate the bead-containing tube by 180°. Wait for the beads to migrate towards the magnet and form a pellet. <input type="checkbox"/> Rotate the tube 180° again (back to the starting position), and wait for the beads to pellet. <p><input type="checkbox"/> Remove the 70% ethanol using a pipette, and discard.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual 70% ethanol.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 20 µl Nuclease-free water. Incubate for 5 minutes at RT.</p> <p><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Pipette 20 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p>In a clean 1.5 ml Eppendorf DNA LoBind tube, mix the reagents in the following order:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 20.0 µl Reverse-transcribed RNA from the "Reverse Transcription" step <input type="checkbox"/> 8.0 µl NEBNext Quick Ligation Reaction Buffer <input type="checkbox"/> 6.0 µl RNA Adapter (RMX) <input type="checkbox"/> 3.0 µl Nuclease-free water <input type="checkbox"/> 3.0 µl T4 DNA Ligase <p><input type="checkbox"/> Mix by pipetting.</p>	

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revK_14Aug2019
 Last update: 27/12/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. <input type="checkbox"/> Resuspend the stock of Agencourt RNAClean XP beads by vortexing. <input type="checkbox"/> Add 40 µl of resuspended RNAClean XP beads to the adapter ligation reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Add 150 µl of the Wash Buffer (WSB) to the beads. Close the tube lid, and resuspend the beads by flicking the tube. Return the tube to the magnetic rack, allow beads to pellet and pipette off the supernatant. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 21 µl Elution Buffer by the gently flicking the tube. Incubate for 10 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 21 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Quantify 1 µl of reverse-transcribed and adapted RNA using the Qubit fluorometer DNA HS assay - recovery aim ~200 ng.</p>	
<p>The reverse-transcribed and adapted RNA is now ready for loading into the flow cell.</p>	
<p>Priming and loading the SpotON flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the RNA Running Buffer (RRB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT before placing the tubes on ice as soon as thawing is complete. <input type="checkbox"/> Mix the RNA Running Buffer (RRB), Flush Buffer (FB) and Flush Tether (FLT) tubes thoroughly by vortexing, spin down and return to ice. <input type="checkbox"/> Open the MinION Mk 1B lid and slide the flow cell under the clip. <input type="checkbox"/> Slide the priming port cover clockwise to open the priming port. <input type="checkbox"/> How to prime and load the SpotON Flow Cell 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µls, and make sure that the array of pores is covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores. 	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles (a few µls):</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip 	

Direct RNA sequencing (SQK-RNA002)

Version: DRS_9080_v2_revK_14Aug2019
 Last update: 27/12/2019



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Prepare the flow cell priming mix: add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing. <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes. During this time, prepare the library for loading by following the steps below. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thoroughly mix the contents of the RRB tube by vortexing or pipetting, and spin down briefly. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Take 20 µl of the prepared RNA library and mix it with 17.5 µl of Nuclease-free water. <p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 37.5 µl RRB <input type="checkbox"/> 37.5 µl RNA library in Nuclease-free water <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of sample to the Flow Cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION Mk 1B lid. 	
<p>Ending the experiment</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Wash Kit instructions and store the washed flow cell at 2-8 °C, OR <input type="checkbox"/> Follow the returns procedure by washing out the flow cell ready to send back to Oxford Nanopore. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol. 	